

PRODUCTION OVERVIEW BITTERFELD

In 1998 CFB branched into API by setting up a state-of-the-art production facility in Bitterfeld in the new German state of Sachsen-Anhalt, Germany. Bitterfeld produces Active Pharmaceutical Ingredients and key intermediates according to cGMP guidelines. The plant was WHO GMP approved in 2000 and FDA inspected in 2002.

Please check our various equipment at our Bitterfeld production facility:

Reaction Vessels

Material:	steel, stainless steel, enameled
Volume:	4,000 to 6,000 liter for production plant
Temperature:	-10° to 150°C
Pressure:	up to 10 bar
Distillation receiver added:	3,000 to 5,000 liter

Solid Separation

Material:	stainless steel, enameled
Preparation:	vertical peeling centrifuge, nutsch
Capacity:	200 kg

Dryers

Model:	paddle, tray and double cone
Material:	enameled, steel, stainless steel
Capacity:	100 to 1,000 kg
Temperature:	up to 90°C
Pressure:	vacuum and normal pressure

Warehouse

Warehouse:	store room for packaging material raw material and intermediate storage finished product store room
------------	---